



HRM-CFD Software Release Notes

Date	19 th April 2012
Rev	3.00
File name	HRM-CFD_Install_3p00.exe
Reason for Change	Change to new "Active Installer". No change to GUI software. Rev to 3.00 to keep in line with other products.

SOFTWARE REVISION HISTORY

CFD GUI Software

Date	Rev	Reason
12-Feb-2010	1.00	First release of code
19-Feb-2011	1.01	1. Change to counting method to allow count logging. 2. Reading of USB was indeterminate. Buffer now purged twice with delay between each before reading.
23-Aug-2011	1.02	Minor cosmetic changes to the dialog box. 1. New LOGO 2. Dialog title is SensL not sensL.

FPGA Firmware

Date	Rev	Reason
4-Feb-2009	1.01	First release of code
19-Feb-2010	1.02	Create new counting method for logging counts

File Names: CFD_1p02.pof
CFD_1p02.rpd

Note: FPGA upgrades should only be carried out if absolutely necessary. If an upgrade is deemed necessary, it should be carried out with SensL supervision.